	ENGLISH	DEUTSCH	ESPAÑOL	FRANÇAIS	ITALIANO	TURKISH
	SYMBOLS USED IN LABELING	BEI DER ETIKETTIERUNG VERWENDETE SYMBOLE	SÍMBOLOS USADOS EN EL ETIQUETADO	SYMBOLES UTILISÉS SUR LES ÉTIQUETTES	SIMBOLI UTILIZZATI NELLE ETICHETTE	ETİKETLEMEDE KULLANILAN SEMBOLLER
	Consult instructions for use	Zum Gebrauch die Gebrauchsanweisung heranziehen	Consulte instrucciones para su uso	Consultez le mode d'emploi avant utilisation	Consultare le istruzioni per l'uso	Kullanım için talimatlara bakın
	Do not re-sterilize	Nicht erneut sterilisieren	No volver a esterilizar	Ne pas stériliser une deuxième fois	Non risterilizzare	İmplantı yeniden sterilize etmeyiniz
	Sterilized using irradiation	Durch Bestrahlung sterilisiert	Esterilizado mediante uso de radiación	Stérilisé par irradiation	Sterilizzato tramite irradiazione	Işınlama ile sterilize edilmiştir
	Do not use if package is damaged	Nicht verwenden, wenn die Packung beschädigt ist	No utilice el kit si el envoltorio está dañado	Ne pas utiliser si l'emballage est endommagé	Non utilizzare in caso di confezione danneggiata	Ambalaj hasar görmüşse kullanmayın
	Non-pyrogenic	Nicht pyrogenes	No pirogénico	Non pyrogène	Non-pyrogenic	Pirojen değildir
	Do not re-use	Nicht wiederverwenden	No volver a usar	Ne pas réutiliser	Non riutilizzare	Yeniden kullanmayınız
	Catalogue number	Katalognummer	Número de catálogo	Numéro de Catalogue	Numero di catalogo	Katalog numarası
	Batch code	Chargencode	Código de lote	Code du lot	Codice batch	Parti kodu
	Use-by date	Verfallsdatum	Fecha de uso ideal	Date limite d'utilisation	Data di scadenza	Son kullanma tarihi
	CE mark and identification number of Notified Body	CE-Marke und ID-Nummer der benachrichtigten Institution.	Marca CE (Comunidad Europea) y número de identificación del ente notificado.	Marque CE et numéro d'identification de l'organisme notifié	Marchio CE e numero identificativo dell'organismo notificato.	CE işareti ve Onaylanan Kuruluşun kimlik numarası.
	Authorized representative in the European Community	Bevollmächtigter in der Europäischen Union	Representante autorizado en la Comunidad Europea	Représentant autorisé dans la Communauté européenne	Rappresentante autorizzato nella Comunità Europea	Avrupa Topluluğu'ndaki yetkili temsilci
	Manufacturer	Hersteller	Fabricante	Fabricant	Produttore	Üretici

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DIRECTIONS FOR USE

Reference the *Cartiva Synthetic Cartilage Surgical Implantation Technique Guide* for further information.

Cartiva SCI is implanted through the use of dedicated accompanying instrumentation designed to provide the surgeon and subject with an implant that is well-seated through a press fit implantation. The implantation procedure is similar to that used for osteochondral autograft or allograft transplantation, where a defect area is removed and resurfaced.

Implantation of Cartiva SCI device has been validated for use with surgical instrumentation distributed by Cartiva, Inc. The Cartiva SCI Instrumentation are available in both REUSABLE and DISPOSABLE configurations which have been validated for their intended function and use with a cannulated drill, and are specific to the size of the device being implanted. The following table references optimal dimensions for successful implantation of Cartiva SCI implants slightly proud (~0.5-1.5mm) with the surrounding cartilage:

Table 1 Cartiva SCI Implant Site Specifications

Lesion size	Implant	REUSABLE Drill Part	DISPOSABLE Instrument Set	Hole Diameter	Hole Depth
Up to 8 mm	CAR-08 / CAR-08-US	MTD-08	MTK-08	7.9 mm	8.3 mm
8 mm to 10 mm	CAR-10 / CAR-10-US	MTD-10	MTK-10	9.5 mm	10 mm

NOTE: For the required instructions specific to the use, handling, cleaning and sterilization of REUSABLE INSTRUMENTATION – reference L20-0088 (US) or L20-0274 (other markets) that accompanies the Cartiva SCI Implant.

DISPOSABLE Instrumentation

The DISPOSABLE instrumentation is provided sterile and may only be used once.

Using standard surgical technique, access the affected joint as necessitated by the size and location of the lesion. Care should be taken to avoid nerve damage along the dorso-medial aspect of the joint. Expose the entire joint to gain access to the central metatarsal head. Resect any osteophytes from the proximal phalanx and/or metatarsal head, ensuring adequate dorsal bone stock is preserved for insertion and stability of the implant. Select the disposable instrument set based on the size of the lesion. Confirm the appropriate size implant to be used by using the concave end of the appropriate Placer size on the metatarsal head.

Once the appropriate size is determined, use the concave end of the Placer, ensuring it is centered in the medial/lateral plane, to determine proper placement of the Placement Guide Pin. Remove the protective caps from the Placement Guide Pin and beware of the sharp tip under the metal protective cap. Insert the Placement Guide Pin into the center of the defect ensuring it is securely seated within the defect and is perpendicular to the central aspect of the metatarsal head. The placer should be positioned relatively central but can be slightly asymmetrical so as to address the worst area of arthritic involvement.

Select the appropriate instrument set (8 mm or 10 mm), and remove the protective caps from the included Drill Bit. Beware of the Drill Bit's sharp blades. Utilize the Drill Bit to drill a hole into the subchondral bone to the proper depth (recommendations can be found in the table above). The set and included drill bit should match the selected implant size to achieve a tight fit with the implant. Drilling should be conducted with the Drill Bit perpendicular to the articular cartilage surface with the Drill Bit centered on the repair area. Insert the Drill Bit over the Placement Guide Pin, and advance the drill until the stop reaches the level of the adjacent tissue, to ensure appropriate depth is achieved. Confirm the cavity depth using the depth mark near the flat end of the Placer.

If necessary, remove any cartilage and/or bone debris from the recipient implant site.

Remove implant from packaging using smooth forceps. Moisten the Introducer tube with sterile saline. Place the Cartiva SCI implant into the wide end of the Introducer with the flat end first (curved portion anterior) so that the flat side of the implant will be placed in the bottom of the joint cavity. Markings on the Introducer depict the proper orientation of the implant in the Introducer. Insert the smaller, flat end of the Placer into the wide end of the Introducer. Rest the distal end of the Introducer on a flat, non-shedding, sterile surface and slowly advance the implant to the distal end of the Introducer using the Placer. Place the distal end of the Introducer into the target implant site. Advance the Cartiva SCI implant into the implant site using the Placer. Remove the Introducer and Placer.

Confirm that the final placement of implant is tight in the implant site. The implant should be slightly proud (~1.5 mm) in the implant site.

PACKAGING

All instruments and implants are provided with dedicated packaging that has been qualified and validated for its particular use.

The Cartiva SCI disposable instrument sets are provided pre-packaged and sterile. They are intended for single use only. Before presentation to the operative field, inspect the package to ensure sterility has not been compromised during transportation. The shelf life of the disposable instrument set is five years. The use-before-date is provided on the shelf box and inner tray labels. Re-sterilization of the disposable instruments is strictly prohibited.

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